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PROTECTIVE GLOVES AGAINST RISKS : CHEMICAL

THIS ITEM HAS BEEN DESIGNED IN COMPLIANCE WITH THE STANDARDS:	EN 420:2003 + A1:2009 : Protective Gloves General Requirements. EN ISO 374-1 : 2016 : Gloves for protection against chemical risks and micro-organisms. EN ISO 374-5 : 2016 : Gloves for protection against dangerous chemical products and micro-organisms.
The CE marking on this glove means that it meets the essential requirements provided for in European Directive EEC 89/686 concerning Individual Protective Equipment (IPE): Safety - Comfort - Dexterity - Sturdiness or Regulation 2016/425 concerning Individual Protective Equipment from April 2018	
This IPE model is subject to a EC test carried out by a certified body:	CTC N° 0075- Parc Scientifique Tony Garnier 4, rue Hermann Frenkel - 69367 LYON Cedex - FRANCE

Category III production surveillance is carried out according to articles 11a or 11b or module C2 or D from April 2018

THE FOLLOWING MARKING APPEARS ON THE BACK OF EACH GLOVE:
 Marking in compliance with the European Directive 89/686 (minimum height 5mm) or Regulation 2016/425

Size of hand (mm)	Around the hand (mm)	Length of hand (mm)
6	152	160
7	178	171
8	203	182
9	229	192
10	254	204
11	279	215


Item reference and size

No. of body ensuring follow-up as per articles 11A or 11B or module C2 or D
If 0197: TÜV Rheinland LGA Products GmbH
 /N°0197/Tillystr.2/90431 Nürnberg/Germany

CE Logo + Symbol

RESULTS & PERFORMANCE LEVELS

Use-by date

KEY TO STANDARDS SHOWN

Mechanical Performance:
 Abrasion Resistance (Cycles)
 Cut Resistance (Newtons)
 Tear Resistance (Newtons)
 Perforation Resistance (Newtons)
 Cut Resistance (TDM) (Newtons)
 0 : Level 1 not reached
 X : Not tested/Not claimed
EN ISO 374-1 (Categorie III chemical)
 Letter in reference to the products used
EN ISO 374-5 micro-organisms

EN ISO 374-5:2016  Virus	Protection bacterie/fungus Impermeability and Penetration Resistance to chemical products and micro-organisms: level 2 (min:0 max: 3) Protection virusus/bacterie/fungus Impermeability and Penetration Resistance to chemical products and micro-organisms: level 2 (min:0 max: 3) ISO16604, Procedure B: no detectable transfer (<1 PFU/ml) of the Phi-X174 bacteriophage in the assay titre	  Virus
Permeation and Degradation Resistance: see enclosed chart or ask your supplier.		
EN ISO 374-1: 2016  XYZ	Permeation level	Permeation time (mn)
	Level 1	t< 10min.
	Level 2	t> 30min.
	Level 3	t> 60min.
	Level 4	t> 120min.
	Level 5	t> 240min.
	Level 6	t> 480min.
	Contact with food: If the marking of the glove includes the food contact symbol, the protective gloves comply with Regulation (EC) no. 1935/2004. For further information concerning the application please refer to the food contact declaration of conformity.	
This glove fulfils the requirements of EC Directive EWG 93/42/EWG and EN 455.		
	Article 742 is compliant in accordance with DIN EN ISO 14644-1: ISO Class 5, GMP A/B	

WARNING: This information does not reflect the actual length of protection in the work place, nor the differentiation between compounds and products composed of pure chemicals. Chemical resistance has been tested under laboratory conditions from samples obtained purely from the palm of the hand (with the exception of those instances where a cuff longer or equal to 400mm has also been checked) and solely involves the chemical product being tested. It may be different if used with a compound. It is recommended to check if the gloves are suitable for their intended use, as the conditions in the work place may differ from those in the test in question, depending on temperature, abrasion and degradation. When used, the protective gloves may offer reduced resistance to dangerous chemical products, due to the alteration in their physical properties. Shifts, tears, scrapes or degradation caused by contact with chemical products etc. may considerably reduce their actual service life. With regard to corrosive chemical products, degradation may be the single most important factor to be taken into account when choosing chemical product resistant gloves. Prior to use, it is recommended that the gloves are inspected to ensure that they do not display any defects. Penetration resistance has been tested under laboratory conditions and solely relates to the item tested. Gloves are for single use only. Not tested against viruses.

DO NOT USE: In hot environments with effects comparable to those of a temperature higher than or equal to 50 C. In cold environments with effects comparable to those of a temperature below than or equal to -5 C. We would like to draw the users' attention to the fact that this glove, with ultra pulling resistance, must not be used when there is a risk of entanglement with machines in motion. The materials and ingredients comprising the gloves are not known to contain substances in such concentrations that they would be suspected of having harmful effects on the health and safety of the user under anticipated working conditions. Certain gloves contain NATURAL RUBBER latex capable of causing allergic reactions. Consult a doctor as soon as possible in the event of an allergic reaction. Claims of metacarpal protection do not apply to fingers. Where this glove is supplied with protection to one part of its surface (a face or a palm), then the levels of protection are guaranteed regarding that part only. This IPE only offers protection against the risks indicated in this technical data sheet and only for the levels indicated. Any risk not included in this technical data sheet is not covered. The analysis of residual risks at the workstation and the choice of the suitable IPE (new or cleaned) is the responsibility of the user (EC

Directive 89/656/EEC).
ANY MODIFICATION OF THIS IPE SHALL LEAD TO THE LOSS OF THE GUARANTEE OF ITS LEVELS OF PROTECTION.
CLEANING : The levels of protection cannot be guaranteed, if cleaned.
STORAGE : Keep flat, away from damp / light, without additional weight on top, in its original packaging at a temperature between 5 °C and 25°C. Do not expose to sunlight or ozone.
These products are not intended to be sold or distributed in California, USA.
The EU Declaration of Conformity for each product is available from our website: https://doc.honeywellsafety.com/

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